

Analysis of Elements in Snack Foods: A Closer Look at Pepperoni, Rice Noodles, Frozen Dinners, and Pizza

Elaine Hasty, Macy Harris, Leanne Anderson, Jenny Nelson, and Greg Gilleland

A set of food samples with different percentage compositions of fats, proteins, and carbohydrates were digested in a single microwave digestion batch and analyzed to determine elemental concentrations using inductively coupled plasma–mass spectrometry (ICP-MS) and ICP-optical emission spectroscopy (OES) using the U.S. Food and Drug Administration (FDA) Elemental Analysis Manual (EAM) 4.7 (for ICP-MS) and FDA EAM 4.4 (for ICP-OES). Built-in software tools were used to streamline the analytical workflow, which is especially useful for new or less experienced users running these methods. Processed foods often contain high concentrations of some elements and trace amounts of others, so a wide analytical range is required. The ICP-MS method used to analyze the varied samples in this work was new. To ensure optimum setup of the ICP-MS method, semi-quantitative results were used to assess the solid content of the food digests before the sample preparation dilution was finalized. The same food samples were analyzed using ICP-MS and ICP-OES and the two methods were compared.

Agricultural food samples contain a wide variety of trace and minor elements, which can vary depending on several factors, including the geographical location, soil type, agricultural practices, and crop type. Along with the elements that have been monitored for decades, newly emerging contaminants may also be present in agricultural soil and other chemicals and additives that enter the food chain. Important emerging contaminants include the technology critical elements (TCEs), such as rare earth elements (REEs) that are used in advanced materials, renewable energy, batteries, consumer electronics, and life science applications. The mining, refining, use, and disposal of these elements may cause environmental contamination, which can ultimately lead to food contamination.

Snack foods contain multiple ingredients and are often highly processed, with each ingredient and process representing a potential source of contamination. Although cases are rare, especially if manufacturers adhere to Current Good Manufacturing Practices (cGMP) (1), processed foods can become contaminated at any point during the production, packaging, and distri-

bution process. Typically, the chemicals that are controlled in foodstuffs include organic contaminants, such as pesticide residue, and inorganic contaminants. According to the general standard for food and feed outlined in the Codex Alimentarius (2) published by the Food and Agricultural Organization (FAO) and World Health Organization (WHO), the most concerning elements include arsenic (As), cadmium (Cd), lead (Pb), and mercury (Hg). The maximum levels (MLs) of each of these elements depend on the type of food, as shown in Table I. Codex MLs ensure food does not contain contaminants at levels which could threaten human health. Many countries base legislation on Codex standards and related texts.

To ensure foods are safe and comply with United States regulations and laws, the U.S. Food and Drug Administration (FDA) develops and publishes works that describe analytical methods that manufacturers or importers can use in their laboratories. Examples include the FDA Elemental Analysis Manual (EAM) 4.4 for inductively coupled plasma–optical emission spectroscopy (ICP-OES) and EAM 4.7 for ICP–mass spectrometry (MS) (3,4). Both ICP-OES and ICP-MS are well-established techniques, and



Environmental & Regulatory



Food, Plant & Animal Tissue



Pharmaceuticals



Inorganic Compounds



Oils & Plastics



Metallurgical

What do you need to digest today? **MARS™ 6 can digest it all.**

With industry leading technologies to make digestion easy and complete, it's no wonder why MARS 6 is the world's best-selling microwave digester. Combine this with round the clock service and application support and you will quickly see why CEM customers are customers for life. At CEM, your success is our goal and we continuously improve our systems to meet your current digestion challenges.

New Samples?

Over 100 years of digestion expertise is built into every MARS 6. Digest anything with hundreds of preprogrammed methods that guide you through sample weight, acid type and volume, and tips to achieve the best results.

Too many samples?

Why have we sold over 10,000 sets of MARSXpress vessels? Because they work! With their simple assembly and ability to digest 40 samples in under an hour MARSXpress continues to meet the demands of high throughput labs worldwide.

CEM understands the needs of busy testing labs and is here to support you along the way. Contact us today to see a MARS 6 digest YOUR samples in YOUR lab and learn why 10,000 CEM customers can't be wrong.

New operators?

Watch training videos right on your system! MARS 6 is simple and intuitive but if you need a refresher on vessel assembly or system features you can learn right on your system.

To much work, not enough time?

Work smarter with iLink remote system operation. Start, stop, and monitor multiple MARS 6 runs from the convenience of a tablet. Take iLink back to your desk so that you can work on other tasks all while you monitor your runs, review previous run data, and gain 24/7 access to CEM technical resources.



We Simplify Science

cem.com/MARS6

© 2019 CEM Corporation. All Rights Reserved.




FIGURE 1: Nutritional composition plot for foods and standard reference materials (SRM) analyzed in this study.

because of their multielement and high sample throughput capabilities, both are widely used in the routine measurement of trace elements in a wide range of samples, including foods and beverages (5–11). EAM 4.4 and 4.7 outline how to determine multiple elements in food digests prepared using microwave-assisted acid decomposition. The methods also outline a series of quality control (QC) tests to ensure instrument performance and data accuracy.

To assist busy analysts who run the EAM methods using ICP-OES or ICP-MS, software tools can speed up method development, improve the usability of both techniques, and provide confidence in the results. Many laboratories that analyze unknown samples or handle a wide variety of samples use semi-quantitative analysis to gain a better understanding of the sample before developing the quantitative method. Knowing more about the samples is useful for all analysts, but software tools that reduce some of the uncertainty of using ICP techniques are especially useful for new users.

Acquiring a full ICP-MS or ICP-OES spectrum adds only a few seconds to the quantitative measurement, but the extra data collected can provide valuable ad-

ditional information on the samples. As a technique, ICP-MS can acquire semi-quantitative results for as many as 78 elements in each sample, which provide a more complete overview of the sample composition and confirmation of the quantitative results. The additional information provides certainty and confidence in the results, preventing the need for remeasurements, potentially saving time and resources over the course of the analytical sequence.

Processed foods often contain high concentrations of calcium (Ca), iron (Fe), potassium (K), magnesium (Mg), sodium (Na), phosphorus (P), manganese (Mn), copper (Cu), and zinc (Zn), whereas trace elements such as As, chromium (Cr), Cd, Pb, nickel (Ni), molybdenum (Mo), selenium (Se), and Hg are required analytes in most regulated methods, so a wide analytical range is required.

Materials and Methods

Samples and Sample Preparation

A set of food samples with different percentage compositions of fats, proteins, and carbohydrates were measured in this study, including pepperoni, rice noodles, frozen dinner, and frozen pizza bought from a supermarket in North Carolina (Figure 1).

The food samples were digested as received, except the frozen pizza and frozen dinner, which were homogenized in a blender before sampling. The samples were prepared for analysis according to the digestion procedure outlined in the EAM method using a MARS 6 iWave closed-vessel microwave digestion system (CEM Corporation). After accurately weighing the samples (approximately 0.5 g of food or SRM) into 75 mL PFA MARSXpress vessel liners, 8 mL of nitric acid (HNO₃) and 1 mL of hydrogen peroxide (H₂O₂) were added to the vessel liners. Duplicates of the samples, standard reference materials (SRMs), and spiked samples were then digested in a single batch using the heating program shown in Table II. Each digestion batch can accommodate up to 40 varied food sample matrices, with a single program being used for all sample types. Finally, 0.5 mL concentrated hydrochloric acid (HCl) was added to the digests, followed by deionized water (DIW) to a final weight of 100 g.

Analytical Sequence

Calibration standards for both ICP-MS and ICP-OES methods were prepared in 1% HNO₃ and 0.5% HCl using standard solutions (environmental calibration standard, multielement calibration standard-1, and 1000 µg/mL single calibration standard for Hg, Agilent Technologies). For ICP-MS analysis, most elements were calibrated from 0.1 to 25 ppb, and Cu, Zn, and Mn were calibrated up to 250 ppb, whereas Hg was calibrated from 0.01 to 2.5 ppb. For ICP-OES analysis, standards were prepared from 10 ppb to 10 ppm for all elements, with a 100-ppm standard for Ca, K, Mg, Na, and P.

For ICP-MS analysis, an internal standard (ISTD) solution containing 2 ppm scandium (Sc), germanium (Ge), rhodium (Rh), indium (In), terbium (Tb), lutetium (Lu), and bismuth (Bi) was prepared in a solution consisting of 1% HNO₃, 0.5% HCl, and 10% isopropanol (IPA). The IPA was added to ensure a consistent level of carbon (C) in the digests and calibration standards. This approach alleviated the concern of signal variability

that could occur for some poorly ionized analytes (such as As and Se) because of ionization enhancement in the presence of C. The ISTD solution was added automatically to the sample and standard solution at a ratio of approximately 1:16 using a T-connector fitted to the sample uptake line that directly lead to the nebulizer.

The ICP-OES ISTD consisted of a 2 mg/L solution of yttrium (Y) in 2% HNO₃, which was added in-line at a dilution of ~1:8. The sample peristaltic pump tubing was white-white (1.03 mm i.d.), and the internal standard was orange-green (0.38 mm i.d.). The standard sample introduction system was used, consisting of a glass-concentric nebulizer and a double-pass, glass cyclonic spray chamber. An independent check solution (ICS) and standard blank were analyzed after calibration, and a check standard and standard blank were analyzed every 10 samples thereafter.

The analytical sequences of calibration standards, samples, and QC solutions for both techniques are shown in Figure 2. Each EAM method specified a minimum number of QC samples to be analyzed with each batch of sample. These QC samples include reference materials (RM), fortified analytical portions (FAP), fortified analytical solutions (FAS), method blanks (MBK), fortified method blanks (FMB), and the number of replicates to be included in each analytical run.

Instrumentation

The accuracy of the semi-quantitative analysis of complex sample matrices by ICP-MS is vastly improved if all elements are acquired using helium collision cell mode to control polyatomic interferences. Helium collision cell mode is widely used to remove polyatomic interferences in multi-element quantitative analysis. Because helium is a nonreactive gas, no new interferences are formed in the cell, whereas common matrix-based polyatomic interferences are removed using kinetic energy discrimination (KED). This makes helium mode ideally suited to the quantitative and semi-quantitative analysis of samples with complex and variable major element compositions (12). In this study, a single quadrupole ICP-MS




FIGURE 2: Analytical sequence for (a) ICP-MS and (b) ICP-OES.

Element	Used	Flags	Wavelength	Rating	Concentration	Intensity	Background
Mn	✓		257.610	★★★★	3.13	5117.1	1931.4
			259.372	★ ?	6.88	6886.8	2760.4
			293.931	★★★	3.26	1164.7	4130.9
			293.305	★★			3513.0

Analyte: Mn(259.372)
Confidence: very weak
Interference: Fe(259.373)
Confidence: very strong

Green Check marks indicate the best wavelength for analysis of this element for this sample.

Hover the question mark to reveal the reason for rejecting this wavelength.

FIGURE 3: Information for manganese (Mn) in a beef jerky sample. The ICP-OES IntelliQuant data automatically produces a star ranking system for wavelengths. This information enables the analyst to easily identify the best wavelengths and what wavelengths to avoid.




FIGURE 4: Periodic table heat map view of ICP-MS semiquantitative data acquired for powdered doughnut, showing a high concentration of titanium.

instrument (7850, Agilent Technologies) fitted with an octopole-based collision–reaction cell (Octopole Reaction System CRC) and Ultra High Matrix Introduction (UHMI) aerosol dilution technology was used.

Although polyatomic ions were successfully removed using helium

cell mode, some other spectral overlaps required a different approach. The REEs have relatively low second ionization potentials, so they can form a small percentage of doubly charged ions (M²⁺) in the plasma along with the normal single-charged (M⁺)

TABLE I: Codex MLs for toxic elements in some different food types (not an exhaustive list)

Food Type	Maximum Level (mg/kg)			
	As	Cd	Pb	Hg
Edible fats and oils	0.1		0.08	
Fat spreads and blended spreads	0.1		0.04	
Canned foods (fruit and vegetables)			(0.1)	
Canned beverages				
Natural mineral waters	0.01	0.003	0.01	0.001
Rice, husked	0.35			
Rice, polished	0.2	0.4		
Cereal grains (wheat)	-	0.1 (0.2)	0.2	
Salt, food grade	0.5	0.5	1	0.1
Cured meats				
Vegetables		0.05 to 0.2	0.05 to 0.3	
Fruits	-	-	0.1 to 0.2	
Chocolate (>70% cocoa)		0.8 (0.9)		
Infant formula			0.01	
Milk			0.02	
Fish			0.3	

TABLE II: Microwave digestion parameters

Power (W)	1800
Ramp Time (min)	25
Hold Time (min)	15
Temperature (°C)	200

TABLE III: Total matrix solids data for four food samples (*n* = 2, in duplicate and corrected for dilution) obtained by the TMS function of ICP-MS software

Conc (ppm)	Beef Jerky	Frozen Dinner	Gummy Bears	Pepperoni	Pizza	Rice Noodles
TMS	36433	5553	768.2	20924	7721	2308

TABLE IV: ICP-MS operating conditions. The parameters are set automatically by the software. Higher energy helium mode gives ppt detection limits for selenium, eliminating the need for a separate reaction gas mode.

ICP-MS Parameter	Setting
RF power (W)	1600
Sampling depth (mm)	10
Carrier gas flow (L/min)	0.6
Dilution (UHMI) gas flow (L/min)	0.35
Lens tune	Autotune
Helium cell gas flow (mL/min)	4.3 (10 for Se)
Energy discrimination (V)	5 (7 for Se)

ions. Doubly charged ions appear at half their true mass in the ICP-MS mass spectrum because the quadrupole mass filter separates ions based on their mass-to-charge ratio (*m/z*). Therefore, REEs, such as neodymium (Nd), samarium (Sm), gadolinium (Gd), and dysprosium (Dy), which have isotopes at masses 150, 156, and 160, will form M²⁺ ions at *m/z* 75, 78, and 80, respectively. If the REEs are present at a high enough concentration in a sample, these REE²⁺ ions can interfere with the measurement of the normal analytical isotopes of As and Se. The EAM 4.7 ICP-MS method recommends that analysts measure ¹⁴⁶Nd, ¹⁴⁷Sm, ¹⁵⁵Gd, and ¹⁶³Dy to monitor the potential overlaps on ⁷⁵As and ⁸²Se. When analyzing unknown or new samples types, the REE content can be quickly and simply assessed using the IntelliQuant semi-quantitative analysis function of ICP-MS MassHunter software (version 4.6 and later). If a high concentration of REEs is expected or identified following the semi-quantitative analysis, a method wizard within the software can set up an automated function to correct for the contribution that REE²⁺ ions make to the signals measured for As and Se.

In this study, the semi-quantitative data were calculated from a full mass quick scan spectrum acquired as part of the quantitative analysis, with only two seconds of additional measurement time. The semi-quantitative results provide valuable information about the elemental content of food samples, including:

- The full elemental composition of each sample. The results can be displayed in a table of concentrations or as a periodic table heat map, enabling a quick overview and easy comparison of each sample's composition (Figure 4).
- Identifying the presence of unexpected major or trace elements not included in the quantitative analysis and identifying contamination during sample preparation. For example, the relatively high concentration of Ti identified in the powdered doughnut from the periodic table heat map shown in Figure 4 might



The Value of Knowing

Agilent Atomic Spectroscopy

Knowing the right information at the right time can help you cut through the everyday complexities of your samples, instrument needs, regulations, other demands on your time and resources – and help you handle unexpected challenges.

The most valuable thing an analytical solution can deliver is the confidence of knowing you'll get the answers you need, when you need them. Driven by insights from customers from around the world, Agilent has created an innovative atomic spectroscopy portfolio – spanning instruments, consumables, standards, services, software, and more – to bring that confidence to your lab.

www.agilent.com/chem/atomic

TABLE V: ICP-OES operating conditions

ICP-OES Parameter	Setting
RF power (kW)	1.20
Replicate read time (s)	10
Pump rate (rpm)	12
Nebulizer flow	0.7
Background correction	Fitted

TMS calculation) uses a mass response profile generated from the quantitative calibration standards. The TMS function is a great diagnostic tool in identifying possible causes of internal standard suppression. It is especially useful when dealing with unknown and potentially difficult food sample matrices by helping the analyst to decide if a sample needs to be diluted or a higher HMI setting is needed. The measured TMS levels for the samples analyzed in this study are shown in Table III.

Based on the TMS levels in these food sample types (Table III), the plasma mode of "HMI-4" was selected. HMI applies an additional flow of Ar gas to dilute the sample aerosol by a factor of four times. This "aerosol" dilution method saves the analyst from having to manually dilute samples,

be unexpected. It was likely because of TiO₂ in the white frosting, but the full mass IntelliQuant spectrum can be used to unequivocally confirm the element from its isotope pattern.

- An estimation of the total matrix solids (TMS) level for each sample run, or for a typical sample of a new sample batch.
- The TMS function uses the semi-quantitative data to calculate the approximate solids levels of a sample.

tive data to calculate the approximate solids levels of a sample.

The calculation excludes gas elements, such as argon (Ar), oxygen (O), and nitrogen (N), together with C, P, sulfur (S), and the halides, ensuring a more accurate result. To simplify method setup and ensure the best possible accuracy for uncalibrated elements, the semi-quantitative data (and

TABLE VI: Mean measured concentrations in three food SRMs corrected for dilution. Mean calculated from two separate digests, each measured twice in triplicate.

Element	NIST 1546a Meat Homogenate				NIST 1549 Non-Fat Milk Powder				NIST 2385 Slurried Spinach			
	Cert. Conc. (mg/kg)	Mean Meas. Conc. (mg/kg)	Recovery (%) [*]	QC ^{**}	Certified Conc. (mg/kg)	Mean Meas. Conc. (mg/kg)	Recovery (%) [*]	QC ^{**}	Certified Conc. (mg/kg)	Mean Meas. Conc. (mg/kg)	Recovery (%) [*]	QC ^{**}
Element (mass) ICP-MS data												
⁵² Cr					0.0026	<DL						
⁵⁵ Mn	0.286	0.285	100	Pass	0.26	0.269	103	Pass	3.81	3.65	96	Pass
⁶³ Cu	0.605	0.602	100	Pass	0.7	0.748	107	Pass	0.90R	0.818	(91)	-
⁶⁶ Zn	17.88	19	106	Pass	46.1	42.4	92	Pass	8.37	7.94	95	Pass
⁷⁸ Se	0.281	0.301	107	Pass	0.11	0.118	107	Pass				
⁹⁸ Mo	0.016R	0.019	(119)	-								
¹¹¹ Cd					0.0005	<DL						
²⁰¹ Hg					0.0003	<DL						
Pb ^A					0.019	0.022	116	Pass				
Element (wavelength, nm) ICP-OES data												
Na (589.592)	47R	47	(87)	-	4970	4915	99	Pass	47	42	89	Pass
Mg (279.553)	368	375.4	102	Pass	1200	1121	93	Pass	178	204.4	115	Pass
K (766.491)	3650	3547.7	97	Pass	16900	15875	94	Pass	2490	2631	106	Pass
Ca (396.847)	624	555.9	89	Pass	13000	12932	99	Pass	624	508	81	Pass
P (213.618)	323.7	346.91	107	Pass	10600	10806	102	Pass				

^{*} FDA Elemental Analysis Manual (Section 3.4 Special Calculations) 3.4 Equation 20

^{**} QC criteria FDA EAM 4.7 (80–120%) for NIST certified values

^R Non-certified reference value

^APb results were based on the sum of the signals measured at mass 206, 207, and 208.

<DL, too low to be quantified

which saves time and reduces the chance of introducing errors or contaminants. When HMI is used, all related settings are auto-tuned for the matrix levels of the target sample types. For ease of setup, other instrument operating settings were optimized automatically using the auto-tune function of the instrument control software. Also contributing to the ease of method setup, all quantitative analytes and the quick scan measurement were acquired in helium mode. Instrument operating conditions are listed in Table IV.

For ICP-OES analysis, semi-quantitative data was acquired using a function of the Agilent ICP Expert software known as IntelliQuant Screening. The software automatically makes wavelength recommendations for analytes based on the elemental composition of a sample and any potential spectral interferences that may arise from the sample matrix. An example for a beef jerky sample is given in Figure 3. The data shows a Fe interference on Mn 259.372 nm, accounting for the high concentration of Mn measured using the 259.372 nm line. In the software (Figure 3), Mn 257.610 nm was given the highest confidence ranking (five stars) so was selected for the quantitative method. Having access to this type of wavelength information makes it easier for the analyst to have confidence in the results that they report.

The ICP-OES sample analysis was conducted per FDA EAM-4.4 for 22 elements using an Agilent 5900 ICP-OES operating in synchronous vertical dual-view (SVDV) mode. SVDV mode allows for both the axial and radial light to reach the detector at the same time. Eliminating the need for separate plasma views shortens analysis times while maintaining low limits of detection (LOD) from axial view and the greater linear range for major cation elements from the radial view. The ICP-OES plasma parameters are shown in Table V.

Results

Quality Control Criteria in EAM 4.4 and 4.7

To verify the sample digestion process, two sets of the three NIST SRMs were analyzed in duplicate using ICP-MS and ICP-OES. As shown in Table VI, the mean concentrations were in agreement with the certified concentrations, meeting the QC criteria requirements of the FDA EAM method of 80–120%. Because not all SRMs are certified for all analytes, blank cells indicate the absence of a certified or reference value. There are no certified values for As, Ni, and thallium (Tl), which are listed in EAM 4.7.

The detection limits in Table VII were calculated from three times the standard deviation of 10 measurements of the calibration blank. Per the EAM methods, minimum QC requirements were included in each analytical run. For ICP-OES, all FAP samples were spiked at 0.05 ppm, and the recoveries are shown in Table VII. For samples with higher naturally occurring concentrations, the 0.05 ppm spike results were not reported. For ICP-MS, the same samples were spiked with all elements at 1 or 50 ppb and measured using ICP-MS. For samples that had naturally occurring elemental concentrations below 5 ppb, a 1 ppb spike is reported. For samples with higher naturally occurring concentrations, the 50 ppb spike results are reported. The recoveries for all elements in the fortified



FIGURE 5: Periodic table heat map view of ICP-OES semi-quantitative data showing relative concentrations of elements (ppm) in beef jerky.

food samples were within the EAM method QC criteria of $\pm 20\%$, as shown in Table VII.

Quantitative Results for Food Samples

Dilution corrected quantitative results obtained using ICP-MS are given in Table VII (beef jerky and gummy bears), and Table VIII shows the dilution corrected quantitative results for frozen dinners, pepperoni, pizza, and rice noodles. In addition to the 12 elements

TABLE VII: Detection limits and concentration and spike recovery results for two food samples. The mean recovery results are based on the analysis of replicate sample digests and each run-in duplicate on the ICP-OES and ICP-MS. *n* = 2.

Element and wavelength (nm)	Detection Limit	Quantitative Results		Fortified Method Blank	Fortified Analytical Portion (Beef Jerky)	Fortified Analytical Portion (Gummy Bears)
		Beef Jerky	Gummy Bears	0.05 ppm spike		
		mg/kg		Mean Recovery ± 1σ (%)		
Al 396.152	0.0009	4.1 ± 0.2	1.4 ± 0.4	104	109	98
As 188.980	0.0016	<DL	<DL	100	103	93
B 249.772	0.0011	0.7 ± 0.1	<DL	109	107	101
Ba 493.408	0.0001	0.7 ± 0.1	<DL	111	104	102
Ca 315.887	0.0006	321.8 ± 14.7	165.1 ± 4.6	109	*	*
Cd 214.439	0.0002	<DL	<DL	103	98	95
Co 230.786	0.0003	<DL	<DL	103	98	97
Cr 267.716	0.0003	0.1 ± 0.0	0.2 ± 0.0	107	101	99
Cu 327.395	0.0003	1.5 ± 0.1	0.1 ± 0.0	107	105	102
Fe 238.204	0.0002	57.8 ± 2.1	2.1 ± 0.0	105	85	97
K 766.491	0.0272	6121.4 ± 322.9	51.1 ± 4.7	*	*	*
Mg 280.270	0.0001	557.9 ± 19.6	5.7 ± 0.1	109	*	105
Mn 257.610	0.0002	2.9 ± 0.2	<DL	107	100	99
Mo 202.032	0.0005	<DL	<DL	109	104	100
Na 589.592	0.0057	14197.4 ± 1718.4	399.3 ± 21.8	90	*	*
Ni 231.604	0.0006	0.4 ± 0.1	<DL	107	102	98
P 213.618	0.0027	3857.8 ± 151.0	35.2 ± 1.5	*	*	*
Pb 220.353	0.0010	<DL	<DL	107	97	98
Sr 407.771	0.0002	1.8 ± 0.1	0.9 ± 0.0	110	100	100
Tl 190.794	0.0096	<DL	<DL	107	97	98
V 292.401	0.0003	<DL	<DL	107	103	99
Zn 206.200	0.0003	88.4 ± 1.6	<DL	102	91	97
Element		µg/kg		1.0 ppb spike		
²⁷ Al	0.091	3700	1390	101 *	117 *	114 *
⁵¹ V	0.078	<DL	<DL	96	104	113
⁵² Cr	0.021	107	162	101	110	102
⁵⁵ Mn	0.095	2630	31.6	99	90	100
⁵⁶ Fe	0.068	54200	1210	98 *	89 *	111 *
⁵⁹ Co	0.002	1.0	1710	101	99	104
⁶⁰ Ni	0.025	259	<DL	102	101	104
⁶³ Cu	0.572	1250	157	100 *	110 *	94
⁶⁶ Zn	0.146	82400	<DL	94 *	88 *	105
⁷⁵ As	0.008	13.2	<DL	93	98	94
⁷⁸ Se	0.031	476	<DL	98	95	103
⁹⁵ Mo	0.093	47.7	7.60	99	98	97
¹¹¹ Cd	0.003	31.6	6.50	97	102	96

*Samples with higher naturally occurring concentrations. ^APb results were based on the sum of the signals measured at mass 206, 207, and 208.

TABLE VII (CONTINUED): Detection limits and concentration and spike recovery results for two food samples. The mean recovery results are based on the analysis of replicate sample digests and each run-in duplicate on the ICP-OES and ICP-MS. $n = 2$.

Element and wavelength (nm)	Detection Limit	Quantitative Results		Fortified Method Blank	Fortified Analytical Portion (Beef Jerky)	Fortified Analytical Portion (Gummy Bears)
		Beef Jerky	Gummy Bears			
Element		µg/kg		1.0 ppb spike		
¹²¹ Sb	0.006	2630	<DL	90	93	89
¹³⁷ Ba	0.012	577	9.70	96	95	106
²⁰¹ Hg	0.003	<DL	<DL	97	84	82
²⁰⁵ Ti	0.070	<DL	<DL	88	88	100
Pb ^A	0.003	9.50	<DL	97	101	102

*Samples with higher naturally occurring concentrations. ^APb results were based on the sum of the signals measured at mass 206, 207, and 208.

specified in EAM 4.7, data is provided for aluminum (Al), vanadium (V), Fe, cobalt (Co), antimony (Sb), barium (Ba), thorium (Th), and uranium (U). Per the Codex elements of concern, the concentrations for As (4.4 to 89.2 µg/kg), Cd (6.5 to 39.4 µg/kg), Pb (6 to 9.5 µg/kg), and Hg (below detection limit) in all snack foods were below any of the Codex MLs (Table I). The high concentration of Al in the pizza sample, which is likely because of an Al food additive, agrees with the literature values (13).

Periodic Table View of Elements in Food Samples

ICP-MS and ICP-OES methods can both be set up to acquire and process semi-quantitative data in addition to quantitative data acquisition. The ICP-MS quick scan function used in this study acquires full mass-spectrum data for every sample with only 2 s of additional measurement time. Because the data is acquired in helium collision cell mode, analytes are free from common polyatomic ion overlaps, ensuring the quality of the results. For ICP-OES, the detector and comprehensive database of wavelengths allows the software to provide real-time sample insights based on each individual spectrum. In both cases, the results can be displayed in a table or as a color-coded periodic table heat map, showing the approximate concentration of all elements present in a sample.

Figure 4 shows an ICP-MS IntelliQuant heat map for powdered doughnuts. The color intensity indicates the concen-

tration ranges of elements in the sample, with a darker color indicating a higher concentration of that element. The semi-quantitative data provides a complete picture of the elements present in the sample, as data can be reported for elements not included in the calibration standards, which is demonstrated by the high concentration of titanium (Ti), labeled as food additive "TiO₂" on the packaging. The semi-quantitative result for Ti in the doughnut samples was ~90 ppm.

Figure 5 shows an ICP-OES generated periodic table heat map with estimated concentrations of elements in a beef jerky sample. The color-coding used to represent the concentration range for the elements is user adjustable.

Conclusions

Like many analytical laboratories, food testing facilities are under pressure to increase productivity while ensuring data quality, compliance, and audit controls, which can be challenging for staff who may be expected to analyze more samples each day or operate instrumentation that is not familiar to them. It is useful if the instrumentation they use is as intuitive as possible and any reference or regulated methods can be implemented quickly, especially by new users.

In this article, we have shown how software tools available for ICP-MS and ICP-OES can help analysts to simplify method setup and reduce the potential for errors. Software functions, such as IntelliQuant,

can also provide valuable information on sample composition and additional analytes beyond the regulated method requirements. These capabilities are demonstrated for the analysis of various foods that were prepared and analyzed in accordance with the U.S. FDA EAM methods 4.4 and 4.7. None of the snack foods were found to contain potentially harmful contaminant elements above Codex maximum levels of concern.

References

- (1) US Food and Drug Administration, *Current Good Manufacturing Practices (CGMPs) for Food and Dietary Supplements* (FDA, Rockville, MD, 2021). <https://www.fda.gov/food/guidance-regulation-food-and-dietary-supplements/current-good-manufacturing-practices-cgmps-food-and-dietary-supplements>.
- (2) Codex Alimentarius International Food Standards, "General Standard for Contaminants and Toxins in Food and Feed, CSX 193-1995," amended 2019, pp. 45–54 (accessed September 2021). www.fao.org/fao-who-codex-alimentarius/sh-proxy/en/?lnk=1&url=https%253A%252F%252Fworkspace.fao.org%252Fsites%252Fcodex%252Fstandards%252FCXS%2B193-1995%252FCXS_193e.pdf
- (3) W.R. Mindak and S.P. Dolan, *U.S. Food and Drug Administration Elemental Analysis Manual, 4.4 Inductively Coupled Plasma-Atomic Emission Spec-*

TABLE VIII: ICP-MS and supplemental ICP-OES quantitative data for food products

ICP-MS All Conc. (µg/kg)	Frozen Dinner	Pepperoni	Pizza	Rice Noodles
²⁷ Al	825.3 ± 68.2	3250.4 ± 152.3	153992.1 ± 275.4	5234.6 ± 679.0
⁵² Cr	75.6 ± 5.1	679.5 ± 67.3	111.1 ± 6.7	140.4 ± 18.6
⁵⁵ Mn	1607.7 ± 66.2	5347.0 ± 213.7	3721.7 ± 24.7	23037.3 ± 147.2
⁵⁶ Fe	4789.5 ± 66.3	13572.5 ± 93.2	24238.1 ± 710.4	12403.7 ± 195.0
⁵⁹ Co	<DL	<DL	<DL	58.0 ± 20.8
⁶⁰ Ni	54.8 ± 2.7	69.2 ± 5.5	123.0 ± 21.5	252.1 ± 36.4
⁶³ Cu	342.7 ± 11.4	794.2 ± 33.6	795.6 ± 16.2	2050.0 ± 23.9
⁶⁶ Zn	4397.1 ± 80.1	25725.5 ± 311.2	10806.1 ± 105.1	18933.5 ± 96.5
⁷⁵ As	4.4 ± 0.6	5.4 ± 0.2	7.8 ± 0.4	89.2 ± 0.9
⁷⁸ Se*	48.3 ± 1.5	349.6 ± 13.5	99.1 ± 2.0	48.0 ± 1.7
⁹⁵ Mo	68.0 ± 1.2	31.9 ± 1.9	204.7 ± 9.0	338.3 ± 4.5
¹¹¹ Cd	24.0 ± 1.1	23.6 ± 0.5	39.4 ± 1.4	28.1 ± 0.3
¹²¹ Sb	<DL	<DL	<DL	<DL
¹³⁷ Ba	487.6 ± 4.0	127.5 ± 6.2	944.7 ± 11.8	1869.7 ± 39.4
²⁰¹ Hg	<DL	<DL	<DL	<DL
²⁰⁵ Tl	<DL	<DL	<DL	<DL
Pb ^A	<DL	7.0 ± 0.1	<DL	6.3 ± 1.6
²³² Th	<DL	<DL	<DL	<DL
²³⁸ U	<DL	<DL	<DL	<DL
ICP-OES All Conc (mg/kg)	Frozen Dinner	Pepperoni	Pizza	Rice Noodles
B 249.772	<DL	<DL	0.3 ± 0.0	0.2 ± 0.0
Ca 315.887	220.7 ± 21.7	218.0 ± 5.4	960.2 ± 21.2	95.4 ± 2.7
Fe 238.204	5.4 ± 0.1	15.1 ± 0.2	26.9 ± 0.8	13.7 ± 0.2
K 766.491	826.8 ± 7.2	2834.2 ± 168.3	1341.5 ± 7.9	1003.4 ± 3.9
Mg 280.270	125.7 ± 0.3	267.2 ± 12.6	247.0 ± 0.5	951.1 ± 3.0
Na 589.592	2129.8 ± 26.1	7660.7 ± 445.4	2233.7 ± 6.1	14.1 ± 0.5
P 213.618	865.7 ± 6.4	1988.8 ± 59.5	1695.8 ± 5.8	2320.1 ± 8.5
Sr 407.771	0.5 ± 0.0	1.0 ± 0.0	1.6 ± 0.0	0.6 ± 0.0
V 292.401	<DL	<DL	<DL	<DL

^A Pb results were based on the sum of the signals measured at mass 206, 207, and 208.

trometric Determination of Elements in Food Using Microwave Assisted Digestion, Version 1.1 (FDA, Rockville, MD, 2010).

(4) P.J. Gray, W.R. Mindak, and J. Cheng, *U.S. FDA Elemental Analysis Manual, 4.7 Inductively Coupled Plasma-Mass Spectrometric Determination of Arsenic, Cadmium, Chromium, Lead, Mercury, and Other Elements in Food Using Microwave Assisted Digestion, Version 1.2* (FDA, Rockville, MD, 2020).

(5) M. Patriarca et al, *J. Anal. At. Spectrom.* **36**, 452–511 (2021).

(6) D. Link, P.J. Walter, and H. Kingston, *Enviro. Sci. Tech.* **32**, 3628–3632 (1998).

(7) R.A. Nadkarni, *Anal. Chem.* **56**, 2233–2237 (1984).

(8) S.P. Dolan and S.G. Capar, *J. Food Compos. Anal.* **15**, 593–615 (2002).

(9) U.S. Food and Drug Administration, *FDA Total Diet Study, March 2021*. <https://www.fda.gov/food/science-research-food/total-diet-study> (accessed October 2021)

(10) S.A. Baker, D.K. Bradshaw, and N.J. Miller-Ihli, *At. Spectrosc.* **20**, 167–173 (1999).

(11) P.J. Gray and W. Cunningham, *J. AOAC Int.* **102**(2), 590–604 (2019).

(12) E. McCurdy and G. Woods, *J. Anal. At. Spectrom.* **19**, 607–615 (2004).

(13) S.M. Saiyed and R.A. Yokel, *J. Food Additives & Contaminants* **22**, 234–244 (2005).

Elaine Hasty, Macy Harris, and Leanne Anderson are with CEM Corporation. **Jenny Nelson and Greg Gilleland** are with Agilent Technologies. Direct correspondence to: jenny.nelson@agilent.com •