

## Application Note for Speciation Analysis

**Species :** Monobutyltin, dibutyltin and tributyltin (MBT, DBT and TBT)

**Sample Type :** Mussel tissue (CRM-477)

**Sample Category :** Biological Tissue

**Analysis :** Species-Specific-Isotope Dilution Analysis+ GC-ICP-MS or GC-MS

### Sample Preparation Procedure

### Recommended Equipment

**Sample Weight :** 0.1g

**Isotopic Spike<sup>1</sup> :** <sup>119</sup>Sn enriched butyltinmix

**Reagent(s) :** 4mL HAcO/MeOH (3:1)

#### Microwave Heating Program

1

**Temperature Target :** 70°C

**Ramp Time :** 1 minute

**Heating Time :** 4 minutes

#### Derivatization Procedure

2

**0,4mL Extract + 4mL buffer HAcO/AcO 0,1 M at pH 4,9**

**Adjust pH to 4.9<sup>2</sup>**

**Add 0,5mL Hexane**

**Add 0,5mL NaBEt<sub>4</sub> (2%w/v)**

**5 minutes mechanical shaking**

#### Clean Up Procedure

3

**Separate Organic Layer**

**Pass Through Florisil Column**

**Elute with 4 mL hexane and collect**

**Preconcentration with N<sub>2</sub> or Ar**

**Florisil Column**

**N-evap or equivalent**

#### Chromatographic Separation

4

**T Injector - 250 °C / Splitless 0.5 min**


**60 °C(0.5min) to 280 °C (4min)**

**Ramp Rate : 30 °C / min**

**Tr-5 Column**

**5% Phenyl Methyl Siloxane**

**30m x 0.25mm i.d. x 0.25 µm coating**



Explorer 24



Explorer 'S' Class

**File: mussel1.doc**

## Detection and Analysis




GC-MS: JEOL JMS-Q1000GC

## Equipment Used



GC-ICP-MS: Focus GC & Thermo X series

## Optimization of Sample Preparation<sup>3</sup>



**Conditions: 70 °C / 4 minutes**

MS MBT  $1071 \pm 79$  ng Sn/g

DBT  $844 \pm 22$  ng Sn/g

TBT  $823 \pm 9$  ng Sn/g

ICP MBT  $1030 \pm 9$  ng Sn/g

DBT  $792 \pm 6$  ng Sn/g

TBT  $841 \pm 9$  ng Sn/g

**Certified Value :**

**MBT  $1013 \pm 189$  ng Sn/g**

**DBT  $785 \pm 61$  ng Sn/g**

**TBT  $900 \pm 78$  ng Sn/g**

<sup>1</sup> Commercially available through ISC-Science (Spain) the addition of the <sup>119</sup>Sn-enriched species to the sample was done in a range of 1-10 times the expected concentration of the endogenous butyltin species in the tissue.

<sup>2</sup> Adjust pH with HCl conc. NH<sub>4</sub>OH conc.

<sup>3</sup> Be aware that the use of species-specific isotope dilution analysis prevents the degradation of TBT during the extraction.