

Application Note for Speciation Analysis

Species : Methyl Mercury

Sample Type : Tuna Fish (BCR-464)

Sample Category : Biological Tissue

Analysis : Species-Specific-Isotope Dilution Analysis + SPME-GC-MS

Sample Preparation Procedure

Recommended Equipment

Sample Weight : 0.1g

Isotopic Spike : ^{201}Hg enriched MeHg¹

Reagent(s) : 4mL TMAH (25%)

1

Microwave Heating Program

Temperature Target : 70°C

Ramp Time : 1 minute

Heating Time : 4 minutes



Explorer 24



Explorer 'S' Class

2

Derivatization Procedure

1mL Extract + 3mL buffer HAc/AcO 0,1M pH 3,9

Adjust pH to 4.0²

Add 1mL NaBEt₄ (2%w/v)

Close SPME vial

Magnetic stirring for 5 min

3

Clean up Procedure

Insert SPME fiber

Sampling

DVD/CAR/PDMS 50µm/30µm

Room temperature during 15 min

4

Chromatographic Separation

T Injector - 250 °C / Splitless 1 min

40 °C(1min), 15 °C/min to 90 °C

30 °C/min to 280°C

Carrier gas: He (1 mL min⁻¹)

DB-5MS Column

5% Phenyl Methyl Siloxane

30m x 0.25mm i.d. x 0.25 µm coating

File: tuna1.doc

Detection and Analysis

Equipment Used



GC-MS: JEOL JMS-Q1000GC

MeHg concentrations obtained in BCR-464 (tuna fish): 4967 ± 87 ng Hg/g³
Certified Value : 5117 ± 160 ng Hg/g

¹In-house synthesized isotopically enriched MeHg: *Appl. Organomet. Chem.*, 2002, 16, 610-615. Addition of the enriched species was done in a range of 1-5 times the expected concentration in the tissue.

²With HCl conc. or NH₄OH conc.

³Be aware that the use of species-specific isotope dilution analysis prevents the degradation of methylmercury during the extraction.