

NEW



Disposable Glass Vessels for Microwave Digestion

Never wash vessels again.

Other digestion systems using disposable vessels are slow, require dangerous pressures, and can only digest half as many samples. The new CEM glass inserts for MARSPress Plus revolutionize disposability for a fraction of the cost.

Increase your lab efficiency and lower your cost per test.


- Weigh sample directly into the disposable insert without static interference
- Up to 24 samples can be digested simultaneously
- Prepare a second set while another batch is running
- Dilute samples directly in the glass insert
- **Discard insert — No washing**



We Simplify Science

cem.com/glass-vessels

Innovative Design



Get clean digests and eliminate vessel cleaning.

Disposable glass vessels are clean enough to achieve <5 ppb recoveries for most elements when run with acid only. Low detection capability may be inhibited due to inherent material contaminants including: B, Na, Mg, Al, K, and Ca.

For more information or to order visit cem.com/glass-vessels or contact info@cem.com.